

Rapid and Accurate Determination of Trace Elements in Copper Concentrate Fused Beads

SolidSample ICPMS – Fully Automated Laser Ablation ICPMS



Matrix Matched LA-ICPMS Calibration Standards:

Linearity and accuracy previously only thought possible through dissolution.

Full Automation:

Continuous robotic sample feed with barcoded sample tracking, data reduction, and result output.

Faster Turn-Around Time:


Weeks of sample prep and analysis reduced to minutes!



Automating Multi-elemental Analysis of Solid Samples by ICPMS using SolidSample ICPMS

- Homogeneous, matrix-matched calibration standards can be used repeatedly for hundreds of analyses before refurbishing for continued use
- Accurate results and linear calibrations: Calibration standards produce R^2 values of ≥ 0.999 and recoveries within $\pm 5\%$ in certified reference materials
- No acid dissolution, acid disposal, or exposure to hazardous materials. No more HF!
- Capable of analyzing >1000 samples in a single experiment
- Automatic quantification and result output powered by iolite data reduction software
- Sample to sample analysis time in <1 minute

SolidSample ICPMS Analytical Cycle



Robotic Arm for Sample Handling

Robotic arm uses vacuum suction cup to deliver samples to, and retrieve samples from rotary carousel.




Rotary Carousel

Carousel rotates to present samples to the analysis position.



SelfSeal Ablation Chamber

Samples are driven up to the self-sealing ablation chamber via piston, creating an airtight seal for solid sample analysis.



Sample and Standard Preparation

A series of Cu Concentrate matrix-matched, Lithium Tetraborate fused beads were created at ESI for calibration standards, using a proprietary method.


Tungsten was added to the beads for internal standardization.

8 spiked calibration standards (S1-S7) were created in accordance with table 1.

OREAS166 Cu Concentrate reference material was fused and treated as an unknown.

Table 1. ESI in-house fused bead calibration standard concentrations.

Calibration Standard ID	Fused Bead Concentrations						
	S1	S2	S3	S4	S5	S6	S7
Concentration (ppm)	0.39	2.3	4.6	9	19	33	50



Experimental Design

The following experiment was designed to:

- 1.) Verify the quality of calibration with ESI in-house Cu Concentrate calibration standards.
- 2.) Demonstrate long term accuracy over a 75 hour run, emulating high-throughput workflow of industrial mining analysis.

Unknown & calibration standards ran every 45 samples

- 1 set of calibration beads (S1-S7) were used repeatedly over the entire experiment (n=20)
- 1 Unknown bead (OREAS166) was used repeatedly over the entire experiment (n=20)

Blank fused beads used to fill entire sample rack.

1000 sample acquisitions

Experimental Conditions

Each sample was treated with a pre-ablation pass to remove any surface contamination.

Each time a bead was sampled, it was ablated with a single line scan.

He cell gas was used in the ICPMS to suppress polyatomic interferences.

Data processing was performed automatically using built-in SolidSample ICPMS data reduction application powered by iolite.

Table 3. SolidSample ICPMS and ICPMS Parameters

SolidSample ICPMS	
Preablation	
Spot size	150 µm
Repetition rate	10 Hz
Fluence	3.5 J/cm ²
Stage speed	150 µm /s
Ablation	
Spot size	100 µm
Repetition rate	40 Hz
Fluence	5.5 J/cm ²
Stage speed	100 µm /s
Transfer gas (He) flow	800 mL/min
ICPMS	
Masses monitored	²⁴ Mg, ²⁷ Al, ⁴⁴ Ca, ⁴⁷ Ti, ⁵¹ V, ⁵⁵ Mn, ⁵⁹ Co, ⁶⁰ Ni, ⁷² Ge, ⁷⁵ As, ⁹⁵ Mo, ¹⁰⁷ Ag, ¹¹⁵ In, ¹¹⁸ Sn, ¹²¹ Sb, ¹³⁷ Ba, ²⁰⁸ Pb, ²⁰⁹ Bi, ²³⁸ U
Dwell time per amu	20 ms
Neb gas (Ar) flow	800 mL/min
Cell gas (He) flow	2.5 mL/min

Results

The resulting data indicates that analytical goals were met. Calibration curves for most elements showed R² values greater than 0.999. Most elements recovered within ± 5% of expected values. Crucial elements for Cu Concentrate analysis, such as As, Ag, Sb, Bi, and Pb showed exceptional results.

Results: Recovery and LOD




Figure 4. Recoveries. Error bars indicate standard error of 20 replicate ablations.

Table 6. Measured vs Accepted Values

Elements	Measured	Long-term Accuracy (ppm)			R2
		Accepted	% Diff	% Recov	
²⁵ Mg	743.6	741.7	0.2	100	0.9997
²⁷ Al	512.4	534.5	-4.2	96	0.9995
⁴⁷ Ti	17.1	18	-5.1	95	0.9989
⁵¹ V	0.4	0.4	1.3	101	0.9999
⁵² Cr	2.7	2.5	7.3	107	0.9997
⁵⁵ Mn	15.6	15.5	0.6	101	0.9993
⁵⁹ Co	151.6	144.1	5.2	105	0.9994
⁶⁰ Ni	4.9	5	-1.1	99	0.9990
⁷² Ga	0.1	0.1	-5.3	95	0.9999
⁷⁵ As	217.3	221.7	-2	98	0.9999
⁹⁵ Mo	0.6	0.6	3	103	1.0000
¹⁰⁷ Ag	0.7	0.7	1.2	101	0.9998
¹¹⁵ In	0.3	0.3	-2.6	97	0.9990
¹¹⁸ Sn	0.4	0.4	4.1	104	0.9995
¹²¹ Sb	1.6	1.6	-1.6	98	0.9999
¹³⁷ Ba	1.6	1.7	-3.1	97	0.9951
²⁰⁸ Pb	10.3	10	3.4	103	0.9999
²⁰⁹ Bi	2.6	2.5	2.5	102	0.9999
²³⁸ U	0.1	0.1	-3.2	97	0.9976

Table 7. Limits of Detection

LOD (Pettke) (ppb)																		
²⁴ Mg	²⁷ Al	⁴⁷ Ti	⁵¹ V	⁵² Cr	⁵⁵ Mn	⁵⁹ Co	⁶⁰ Ni	⁷¹ Ga	⁷⁵ As	⁹⁵ Mo	¹⁰⁷ Ag	¹¹⁵ In	¹¹⁸ Sn	¹²¹ Sb	¹³⁷ Ba	²⁰⁸ Pb	²⁰⁹ Bi	²³⁸ U
727.9	496.0	55.1	4.8	31.9	33.4	5.6	233.8	4.0	23.4	3.7	4.2	5.0	332.1	2.9	3.8	1.5	0.3	0.3

Results: Calibration Curves






Figure 5a. Calibration Curves

Results: Calibration Curves


⁵⁹Co




⁶⁰Ni




⁷¹Ga



⁷⁵As



⁹⁵Mo



¹⁰⁷Ag




Figure 5b. Calibration Curves

Results: Calibration Curves










Figure 5c. Calibration Curves


Conclusion

SolidSample ICPMS has been demonstrated to provide accurate and precise data for trace elements Cu Concentrate fused bead samples.

ESI Cu Concentrate calibration standards accurately quantified OREAS166 Cu Concentrate samples.

- Calibration curves for most elements showed R^2 values of better than 0.999.
- Almost all elements recovered within $\pm 5\%$ of known values.

SolidSample ICPMS can meet the rigorous demands for high throughput mining samples with a lower cost of ownership compared to other techniques.



Robotic Sample Handling

- Fully automated sample handling
- Unattended 24/7 operation
- High sample transport efficiency
- Process tens to thousands of samples per day
- Up to 5X faster than conventional laser ablation

Laboratory Safe

- Solid sampling eliminates hazards associated with acid digestion
- Class 1 laser-safe during operation

Advanced ICPMS Interface

- Ensures plasma stability during rapid sample exchanges
- Provides exceptional transport efficiency

Sample Carousel

- Swift sample changeover
- Delivers samples to ablation chamber for analysis

Barcode Scanner

- Ensures data integrity
- Tracks sample through the complete analysis workflow

Figure 1. SolidSample ICPMS features diagram.



© Elemental Scientific Lasers LLC | 685 Old Buffalo Trail | Bozeman, MT 59715
Tel: 406-586-3159 | lasers@icpms.com | www.icpmslasers.com

UA-25099