



SolidSample ICPMS

Robotic Sample Handling

- Fully automated sample handling
- Unattended 24/7 operation
- High sample transport efficiency
- Process tens to thousands of samples per day
- Up to 5X faster than conventional laser ablation

Laboratory Safe

- Solid sampling eliminates hazards associated with acid digestion
- Class 1 laser-safe during operation

Advanced ICPMS Interface

- Ensures plasma stability during rapid sample exchanges
- Provides exceptional transport efficiency

Sample Carousel

- Swift sample changeover
- Delivers samples to ablation chamber for analysis

Barcode Scanner

- Ensures data integrity
- Tracks sample through the complete analysis workflow

